Socket Cone

Art. Nr.:
P527871 (Medium)
P527872 (Large)
P527873 (X-Large)

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Neuburger Str. 35, 90451 Nürnberg
Tel. 0911 / 64339 – 10, Fax 0911 / 6492253

English English English English English

The purpose of use:

The Socket Cone complies of a low temperature thermoplastic, which is thermo plastically re-formable through the application of 70°-80°C warm water. It serves for copying and also as a socket material for patients after below knee amputation. Only a qualified specialist (CP / CPO) is allowed to work with the Socket Cone. It is essential that the Socket Cone will be checked after 6 months period of time. The Socket Cone is permitted for the care of only one patient.

Indication:

Below knee amputation.

Application instructions:

1. The choice of Socket Cone

- Measure the circumference of the limb (with donned liner) distal and in a height of 20 cm.

| Size | Circumference (distal) | Circumference (height 20 cm) | Item.-No.: |
|------|---------------------------|---------------------------------|------------|
| M | 18 cm | 22 cm | P 527871 |
| L | 23 cm | 30 cm | P 527872 |
| XL | 28 cm | 38 cm | P 527873 |

2. Preparation of the patient:

- Apply a liner to your patient according to the liner - instructions.
- Screw a pin into the liner's distal attachment.
- Subsequently, isolate the liner with a cling film.

3. Preparation of the Socket Cone:

- Mark the required length (stump length) and warm up this area in a 70°-80°C warm water bath.
- Cut the Socket Cone to the length of your mark with a pair of scissors and apply some Vaseline on the inside and on the outside of the Socket Cone.

- Warm up the whole Socket Cone in a warm water bath with a temperature of 70°-80°C until it becomes completely flexible.
- Pull out the flexible Socket Cone from the warm water and let it cool down slightly.
- Pull the Socket Cone onto the patients limb over the isolated liner until it is evenly covered.
- Attention: The distal drilling Ø 10mm of the Socket Cone must not be stretched.
- Form the Socket Cone with both hands according to the required shape (Hybrid, KBM, PTB etc.) and let it cool down for a few minutes on the patient's leg.
- Now, remove the Socket Cone and dip it in cool water in order to cool it down completely.
- Cut the proximal edge of the socket as desired. Use leather scissors or an oscillation saw or a jig saw and polish it with a sharp grinder operated in low speed.
- Drill the distal hole in the Socket Cone so that it matches to the size of the funnel of your lock system.
- Put your Clutch-Lock incl. dummy for the locking mechanism into the drilling. Make sure that the funnel does not over lap on the socket's inner side.
- Fill all the openings with plasticine and fix a collar made of Trolen on the distal end of the socket.
- Fill the collar with hard foam (H300).
- Sand down a distal area for the alignment of the socket (Flexion / Extension, Abduction / Adduction). Glue your socket attachment (with sealing resin) onto this area in the static as required by your patient.

- After hardening of the glue, grind the foamed area, including the socket attachment according to the socket's contour.

- Go through a dynamic fitting with your patient. Connect your socket with parts chosen for your patient and swap the dummy thread with the locking mechanism and check the function of the pin.

Important: During the dynamic fitting you should make all the necessary amendments to the socket and the build up. After the completion of the socket the possibilities for modification will be limited.

- If the manufacturing of the socket was successful, then it should fit the patient well, allowing some slight form corrections, (e.g. fibular head).

4. Completion of the Socket Cone

- Remove all the modular parts, fit the laminating dummy into the four - hole plate and swap the locking mechanism with the dummy thread.
- Sand down the socket's interior with a rotary rasp head as thin as possible.
- Seal the opening for the pin of your lock system with plasticine.

- Extend and seal the proximal edge of the socket with a plaster loguette.

- Fill the socket with sand, up to about 5cm below the socket's edge (closed area).

- Next, pour plaster into the socket and place a tube in the middle of the socket.


- After the plaster is hardened and cooled, remove the plaster loguette and isolate the whole foamed area including attachments and laminated dummy as well as the proximal plaster seal with sealing resin.

- Put the isolated model into a vacuum pipe and reinforce this model according to the patient's requirements (body weight and mobility level).

- After completing the reinforcement, cover the model with a PVA-bag, seal it on the vacuum pipe and set it under vacuum.
- Fill in the required quantity of resin (according to the volume of the socket) and then seal the fill in area.
- Laminate the socket in the usual way. Make sure that your reinforcement is saturated well and that there are no air bubbles.
- After the laminate is hardened and cooled, cut the proximal edge of the socket, using an oscillation saw and remove the vacuum pipe.
- Remove the plaster seal and subsequently sand as well as plasticine from the socket.
- Grind the laminated dummy to expose the four holes plate. Stop grinding as soon as you reach the plain surface of the dummy. You will also achieve a sided trimming for the socket adaptor.
- Grind the edge of the socket with a fine but sharp sanding drum operated in low speed. Subsequently achieve a very smooth and clean surface by using a fine sanding paper and sanding by hand.
- Expose the tread dummy by grinding. Stop grinding as soon as you have reached the plain surface of the dummy. By doing so, you will achieve a complete trimming of the locking mechanism with an exception of the push button.
- Remove the thread dummy and mount the locking mechanism.
- At the end, warm up the socket's inner wall very carefully with the use of a heat gun. Through this warming, any harsh lines or creases which possibly occurred by isolating the liner during the forming of the Socket Cone will disappear.

Service:

During the half yearly routine examination, the Socket Cone and its attachments should be tested for any possible damages as well as for its firm fitting.

Liability for Material Defects:

Guarantee applies according to our general terms and conditions of business, which have been handed over to the customer prior to closing the contract of sale. A liability for materials is not effective for damage, caused by improper use, lack or wrong care or forming to wrong fit. The use-by period essentially depends on proper application and care.

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